




RESTORE THE CERVICOMENTAL ANGLE


1. **Simple crossing** - of the threads at the centre of the cervicomental angle.
2. **Superficial implantation** - in the subcutaneous tissue, and very superficial over the 5 cm of the middle section.
3. **Tacking the skin backwards** - to restore the cervicomental angle.
4. **Slight tension** - to avoid a feeling of strangulation.




INFINITE-THREAD® IMPLANTATION LAYERS




under the skin / scalp



under the fascia



in the subcutaneous fat



in the lower part of the skin